



Kuraray CLEARFIL bonding agents are the ideal dental adhesives. Kuraray dental products include dual-cure, light and self-cure, chemical-cure, and light-cure bonding agents and systems.

WHY CLEARFIL?

- Radically reduced post-operative sensitivity & improved delivery methods
- Easier than ever application methods
- Specially formulated to solve specific procedural and case needs
- Proprietary MDP monomer is the technological breakthrough behind this award winning product line
- Virtually no post-operative sensitivity with any Kuraray bonding agents

Clearfil BOND SYSTEMS

UNIVERSAL BOND ~ 1-step system, 10 sec light cure, direct restoration, self-etching, total etching or selective etching procedures, 5-10 micron film thickness

SE BOND 2 ~ 2-step system, 10 sec light cure, direct restoration, 50 micron film thickness

S³ BOND + ~ 1-step system, 10 sec light cure, direct restoration, 5-10 micron film thickness

SE PROTECT ~ 2-step system, 10 sec light cure, direct restoration, 50 micron film thickness

DC BOND ~ 1-step system, dual cure, direct restoration, 5-10 micron film thickness

Photo BOND ~ 2-step etching, dual cure, direct restoration, indirect restoration, combination self-cure composite, 1-micron film thickness

DC Activator ~ Can be used with any of the above bonding agents to take them from light-cured to dual-cured



To Order Call Meg
@ 800-992-8987 x4750
DERMATEC DIRECT
4430 E. Adamo Dr. STE. 306 TAMPA, FL 33605

